ACIST Navvus® II


[en] ENGLISH Navvus® II MicroCatheter Instructions for Use	5	[nl] NEDERLANDS Navvus® II-mikrokatheter Gebruiksaanwijzing	37
[cs] ČEŠTINA Mikrokatétr Navvus® II Návod k použití	9	[no] NORSK Navvus® II-mikrokateter Bruksanvisning	41
[da] DANSK Navvus® II MicroCatheter Brugsanvisning	13	[pl] POLSKI Mikrocewnik Navvus® II Instrukcja użytkowania	45
[de] DEUTSCH Navvus® II-Mikrokatheter Gebrauchsanweisung	17	[pt-BR] PORTUGUÊS DO BRASIL MicroCateter Navvus® II Instruções de uso	49
[es] ESPAÑOL Microcatéter Navvus® II Instrucciones de uso	21	[ru] РУССКИЙ Микрокатетер Navvus® II Инструкции по применению	53
[fi] SUOMI Navvus® II -mikrokatetri Käyttöohjeet	25	[sk] SLOVENSKY Mikrokatéter Navvus® II Pokyny na použitie	58
[fr] FRANÇAIS MicroCathéter Navvus® II Mode d'emploi	29	[sv] SVENSKA Navvus® II mikrokateter Bruksanvisning	62
[it] ITALIANO Microcatetere Navvus® II Istruzioni per l'uso	33	[tr] TÜRKÇE Navvus® II Mikrokateter Kullanım Talimatları	66

[en] Symbols [cs] Symboly [da] Symboler [de] Symbole
 [es] Símbolos [fi] Symbolit [fr] Symboles [it] Simboli
 [nl] Symbolen [no] Symboler [pl] Symbole
 [pt-BR] Símbolos [ru] Обозначения [sk] Symboly
 [sv] Symboler [tr] Semboller

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Instructions for Use

REQUIRED MATERIALS

Each ACIST Navvus® II MicroCatheter package contains the following components:

- one single-use disposable microcatheter
- one dispenser coil
- one monorail stylet

Other materials that are required but not supplied include the following:

- 5F – 8F guide catheter
- Standard (0.014 in (0.36 mm)) guidewire
- hand syringe used for flushing the dispenser coil and catheter
- sterile saline solution

CAUTIONS

- Federal (USA) law restricts this device to sale by or on the order of a physician.
- Carefully read all instructions in this document, as well as in the ACIST RXi® System User's Guide or the ACIST Mini™ User's Guide, prior to use. Observe all contraindications, warnings, and precautions noted in both documents. Failure to do so may result in patient complications. The ACIST Navvus II MicroCatheter should be used by physicians with adequate training in the use of this device. ACIST relies on the physician to determine, assess, and communicate to the patient all foreseeable risks of the procedure.
- The device is provided sterile for single use only. Reuse of a single-use device creates a potential risk of patient or user infection. Contamination of the device may lead to illness or serious patient injury. Do not use the product beyond the expiration date indicated on the labeling.

DEVICE DESCRIPTION

The ACIST Navvus II MicroCatheter is a single-lumen monorail catheter designed to be used with standard 0.014 in (0.36 mm) guidewires in the arterial vasculature. The microcatheter is compatible with the ACIST RXi family of system hardware which includes the RXi System and RXi Mini. Features unique to each system are denoted specifically.

The Navvus II MicroCatheter distal shaft is 26 cm in length with a pressure sensor located 5 mm from the distal tip. The elliptical distal shaft is 1.68 x 1.91F (0.020 in x 0.025 in) up to 10 mm from the distal tip; a maximum profile of 2.7F (0.036 in) occurs at the pressure sensor. The distal shaft smoothly tapers over the pressure sensor and down to the tip accepting the guidewire. A radiopaque marker band is located 2.5 mm from the distal tip.

The shaft proximal to the monorail section is 2.4F, allowing use in 5F or larger guide catheters. Two white positioning markers are located at 80 cm and 100 cm from the distal tip.

The pressure sensor on the catheter utilizes optical sensing technology. Both the optical pressure signal and information for auto calibration are transmitted from the catheter to the RXi hardware.

ACIST RXi System: The sensor pressure signal is processed by the console, which displays values in real time.

ACIST RXi Mini: The sensor pressure signal measured by the Navvus Catheter is sent from the processing unit of the RXi Mini to the hospital's hemodynamic system in real time.

The ACIST Navvus II MicroCatheter is not made with natural rubber latex.

INTENDED USE / INDICATIONS FOR USE

The ACIST RXi System and the RXi Mini are indicated for obtaining intravascular pressure measurements for use in the diagnosis and treatment of coronary and peripheral artery disease. The ACIST Navvus II MicroCatheter is intended for use with the entire family of ACIST RXi Systems.

CONTRAINDICATIONS

The ACIST Navvus II MicroCatheter is contraindicated for use in the cerebral vasculature.

INTENDED PATIENT POPULATION

This product is designed for use in adult patients identified by a physician as suitable candidates for vascular pressure measurement, taking into consideration the patient's anatomy and health status.

WARNINGS

The following warnings refer to hazards that can cause serious injury or death.

- Prior to use and whenever possible during the procedure, carefully inspect the Navvus II MicroCatheter for kinks or any other damage. Do not use a kinked or damaged microcatheter as vessel damage and/or the inability to advance or withdraw the microcatheter may occur.
- The product is supplied sterile in a protective carton containing a sterile-barrier tray. Sterility is indicated by a chemical process 'dot' label, which during successful sterilization turns green. Verify package integrity before use. If the sterile barrier package is damaged or the seal is not intact, or if the dot is not green, DO NOT use the product. Contact an ACIST representative to return damaged product. Do not use this product beyond the 'Use by' date indicated, as product integrity and sterility cannot be guaranteed.
- Any serious incidents related to the product or its use must be reported to ACIST Medical Systems, Inc., as well as the competent authority of the Member State or country in which the user and/or patient is established.
- Improper technique may damage the catheter. Exercise care in handling the Navvus II MicroCatheter during a procedure to reduce the possibility of accidental breakage or kinking.
- When delivering the Navvus II MicroCatheter over the guidewire, ensure that the guidewire and the microcatheter move freely within the vessel wall. Failure to do so may traumatize the vessel.
- Use proper technique when crossing lesions.
- The Navvus II MicroCatheter is not designed to be torqued. Do not excessively torque the microcatheter.
- Never advance or withdraw the Navvus II MicroCatheter against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the microcatheter or guidewire against resistance may result in separation of the microcatheter or guidewire tip, damage to the microcatheter, or vessel perforation.
- Once inserted, do not remove the purple connector from the system during the procedure. If using the RXi System, removing it ends the case. If using the RXi Mini, the pressure signal stops transmitting until the connector is plugged back in.
- Store the Navvus II MicroCatheter's working length in saline anticoagulant solution between uses with the same patient.
- This product should not be used in rooms containing magnetic resonance imaging (MRI) equipment.

CAUTIONS

The following precautions refer to hazards that could result in minor injury or damage to the system or other equipment.

- To avoid inaccurate aortic pressure measurements, the use of guide catheters larger than 8F or guide catheters with side holes are not recommended. Avoid blocking blood flow to the selected artery. Do not wedge the guide catheter in the selected artery.
- The Navvus II MicroCatheter is not compatible with 4F guide catheters.
- Carefully place the Navvus II MicroCatheter on the guidewire to avoid damage to the tip. Avoid contact between the guidewire and the small pressure sensor just proximal to the tip.
- Open the hemostasis valve fully to avoid damage to the Navvus II MicroCatheter tip.
- Always advance or withdraw the Navvus II MicroCatheter slowly and carefully.

- Do not perform high pressure (greater than 600 psi) fluid injections while the Navvus II MicroCatheter tip is inside the guide catheter.
- Observe all Navvus II MicroCatheter movement in the vessel. When the microcatheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the microcatheter without observing the resultant tip response.
- Tighten the Tuohy Borst adapter on the guide catheter after placing the Navvus II MicroCatheter distal shaft inside the guide catheter, or false low aortic pressure readings may result.
- Do not move the Navvus II MicroCatheter or the aortic transducer during equalization or recording.
- Confirm Pd/Pa=1.00 after equalization before moving the Navvus II MicroCatheter from the equalization site. After a measurement, retract the microcatheter tip to the equalization position and verify Pd/Pa=1.00.
- Ensure proper placement of the sensor with respect to the lesion being evaluated. Do not start recording until the Navvus II MicroCatheter is properly positioned.
- Once connected, do not disconnect the Navvus II MicroCatheter from the system during the procedure.

PACKAGING AND STORAGE

The Navvus II MicroCatheter has been sterilized with ethylene oxide. Store the Navvus II MicroCatheter at ambient conditions.

POTENTIAL ADVERSE EVENTS

Potential complications that may be encountered during all catheterization procedures include, but are not limited to: vessel dissection or occlusion, perforation, embolus, spasm, local and/or systemic infection, intimal disruption, distal embolization of blood clots and plaque, myocardial infarction, serious arrhythmias, or death.

INSTRUCTIONS FOR USE



The following instructions provide technical direction, but do not eliminate the need for formal training in the use of the Navvus II MicroCatheter and the ACIST RXi system. The techniques and instructions provided here do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.

Before unpacking and preparing the Navvus II MicroCatheter, ensure that the RXi System is on.

1. Carefully inspect the Navvus II MicroCatheter packaging and components for damage prior to use.
2. Using sterile technique, open the Navvus II MicroCatheter package as shown:




Technician:

- a. Place your left thumb on the plastic corner just under the "Open" symbol and grip the underside of the tray at the corner with your fingers.
- b. Hold the tray securely with your left hand and use your right hand to peel the cover back and off the tray.
- c. Present the opened tray to the scrubbed and sterile user. Continue to hold the opened tray for the scrubbed and sterile physician until the product is removed from the tray as described below.


Physician/Technician:


- d. Grasp the purple handle assembly.
- e. Using the same hand and while still holding the purple handle, grasp the orange coil.
- f. Using the same hand and while still holding the purple handle and orange coil, finally grasp the dispenser coil.

3. Place the dispenser coil and the handle assembly cable in the sterile field. Ensure the handle assembly does not get wet.
4. Attach the catheter's handle assembly to the ACIST RXi system as shown below:

If using the ACIST RXi Mini, attach the catheter's handle assembly as shown below:


Notes:

- When manipulating the Navvus II MicroCatheter, avoid creating a knot in the cable.
- The pressure sensor on the Navvus II MicroCatheter is automatically zeroed when the handle assembly is inserted into any RXi interface. Therefore, the tip should be outside the patient when this connector assembly handle is inserted into any RXi hardware.
- The Navvus II MicroCatheter status indicator LED illuminates GREEN when the handle assembly is properly inserted. If the LED illuminates RED, re-insert the handle assembly or replace the microcatheter.

- Using sterile technique, attach a new, sterile hand syringe filled with sterile saline solution to the flushing Luer entry port of the dispenser coil.

Thoroughly flush the Navvus II MicroCatheter, as shown:

- Using sterile technique, remove the Navvus II MicroCatheter from the dispenser coil and inspect it for any kinks. If damage or kinks exist, discard the microcatheter and use another.

To position the ACIST Navvus II Catheter:

The following steps assume a standard PCI protocol using the following items: an inserted 5F – 8F guiding catheter, an inserted 0.014 in (0.36 mm) guidewire, and a hemostasis valve.

Use proper anticoagulation and anti-platelet therapy before beginning.

Note:

Ensure that the RXi hardware is on and cables are connected correctly.

Carefully read all warnings and cautions before beginning.

- Remove the stylet from the end of the Navvus II MicroCatheter. **Do not discard.**
- Advance the Navvus II MicroCatheter over the guidewire that is already in place in the distal vasculature.

Note:

- Carefully place the Navvus II MicroCatheter on the guidewire to avoid damage to the tip. Avoid contact between the guidewire and the small pressure sensor that is just proximal to the tip.
- Open the hemostasis valve fully to avoid damage to the Navvus II MicroCatheter tip.

- On the guide catheter, tighten the Tuohy Borst adapter. Failure to do so may result in a false low aortic pressure reading.
- Advance the Navvus II MicroCatheter over the wire and position it with the radiopaque marker in the coronary ostium, 3 mm outside the tip of the guide catheter.

- Equalization must take place with the Navvus II MicroCatheter distal tip just distal to the guide catheter tip. Pressure is measured less than 3 mm proximal from the marker band visible under fluoroscopy.
- Flush the guide catheter with saline solution.
- ACIST RXi System:**
 - Verify that the Pa displayed on the RXi console matches the reading on the hemodynamic monitor.
 - Press the **Equalize** button.
 - After equalization, the red and green waveforms should be superimposed and the displayed Pd/Pa should be stable at 1.00.
- ACIST RXi Mini:**
 - Equalize the Aortic and Distal pressures with the hemodynamic system FFR software.
 - After equalization, the Aortic and Distal waveforms on the hemodynamic system should be superimposed and the displayed Pd/Pa should be stable at 1.00.
- Advance the Navvus II MicroCatheter sensor past the lesion under study.

PRESSURE PERFORMANCE SPECIFICATION

Operating pressure	-30 to 300 mmHg
Zero drift	≤ 7 mmHg over one hour
Total accuracy	3 mmHg or 3% of reading

ACIST Navvus® II MicroCatheter Limited Warranty

ACIST Medical Systems, Inc. (“ACIST”) warrants that reasonable care has been used in the design and manufacturing of the NAVVUS II intravascular pressure sensing microcatheter (the “Product”). However, handling, storage, cleaning, and sterilization of the Product, as well as factors relating to the patient, including but not limited to diagnosis, treatment, surgical procedures, and other matters beyond the control of ACIST, may affect the Product and/or the results obtained from its use. Accordingly, the obligation of ACIST under this warranty is limited to the repair or replacement of the Product, and ACIST shall not be liable for any direct, incidental or consequential loss, damage, or expense directly or indirectly arising from the use of the Product. ACIST neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with the Product.

The Product is intended for use by a physician. It is important to read the accompanying instructions for use, as well as any associated warnings and potential complications, prior to using the Product.

This warranty is in lieu of and excludes all other warranties not set forth herein, whether express or implied by operation of law or otherwise, including but not limited to, any implied warranties of merchantability or fitness for a particular purpose. The Product is intended for single use only, and ACIST assumes no liability with respect to Product that have been reused, reprocessed, or re-sterilized, and makes no warranty whatsoever with respect to such reused, reprocessed, or re-sterilized Product.

ACIST RXi SYSTEM ONLY: STUDY DESIGN INFORMATION

A prospective, observational, multi-center study of the RXi System and Navvus Catheter was conducted. The primary objectives of the study were to assess the relationship in FFR measurements between the Navvus catheter and the St. Jude PressureWire™ System (St. Jude System), and secondly, the rate of device success defined as a valid FFR measurement by the Navvus Catheter where a valid FFR measurement was taken using the St. Jude System. The study used paired measurements from the Navvus Catheter and the St. Jude System within subjects as controls. Pressure measurements were assessed intra-operatively; subjects were evaluated for safety events through hospital discharge.

Subject Enrollment

Three centers participated in the study and enrolled 58 subjects. Patients with single or multi-vessel coronary artery disease indicated for FFR and a vessel reference diameter of ≥ 2.5 mm were eligible.

The average subject was a 65-year-old male (80%), ranging in age from 42 to 86 years. The majority of subjects had a history of dyslipidemia (86%), hypertension (74%), and angina (64%).

Study procedures were performed through the femoral (32%) or radial (68%) artery. Reference vessel diameters ranged between 2.5 and 5.5 mm (average 3.2 ± 0.6 mm) with percent diameter stenosis ranging from 40 – 90% (average $57.2 \pm 10.0\%$).

Results of Clinical Investigation

- The Navvus Catheter was 100% successful in taking valid FFR measurements in all lesions in which the St. Jude System was able to take valid FFR measurements.
- The overall bias, as determined from a Bland-Altman analysis was -0.01
- Using a Bland-Altman analysis, the observed 95% limits of agreement were -0.13 to 0.10.
- The Kolmogorov-Smirnov test for linearity showed a cumulative p-value of 0.72, meeting success criteria of being > 0.1 .
- The intercept determined by a Passing-Bablok linear regression between FFR measurements made using the Navvus Catheter and that made using the St. Jude System showed a 95% confidence interval of -0.15 to 0.19, meeting success criteria of the interval containing 0.
- The slope determined by a Passing-Bablok linear regression between FFR measurements made using the RXi System and that made using the St. Jude System showed a 95% confidence interval of 0.76 to 1.17, meeting success criteria of the interval containing 1.
- Receiver Operator Curves, using the first measured St. Jude FFR measurement as the standard, created for both the Navvus Catheter measured FFR and the last St. Jude FFR demonstrated very similar areas under the curve (87.4% and 88.1%, respectively).
- The RXi System demonstrated a statistically significant lower rate of drift compared to the St. Jude System, both for mean drift (0.02 ± 0.02 for the RXi System vs. 0.06 ± 0.12 for the St. Jude System) and clinically relevant drift (13% vs. 33% rate of occurrence).

This clinical study also demonstrated a very good safety profile. Serious adverse events occurred during five subjects' hospitalizations; however, no SAEs were related to the RXi System or Navvus Catheter. Similarly, none of 17 reported adverse events were device related. Eleven of the 17 were related to the procedure, including hematomas, mild hemorrhage, elevated cardiac enzymes, nausea, vasovagal reaction, headache, hypokalemia, temporary rhythm disturbances, and wheeze. All the adverse events were mild to moderate in nature.

Proprietary Information Notice

This document contains information proprietary to ACIST Medical Systems, Inc. All rights reserved.

Web site (worldwide): acist.com

Patents: acist.com/patents

No part of this document may be reproduced, transmitted, processed, or recorded by any means or form, electronic, mechanical, photographic or otherwise, nor be released to any third party without the written consent of ACIST Medical Systems, Inc.

© 2019 ACIST Medical Systems, Inc. All rights reserved.

The written and graphic product descriptions in this manual were effective at the time of printing. ACIST Medical Systems, Inc. reserves the right to change specifications and designs without prior notification.

ACIST, Navvus, and ACIST RXi are registered trademarks of ACIST Medical Systems, Inc. RXi Mini is a trademark of ACIST Medical Systems, Inc.

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Bracco Pty, Ltd
ACIST Shanghai Trading Co. Ltd.